

Washers & Sanitizers to automatically wash and apply disinfectants to wheelchairs and other hard-to-clean medical equipment.

# AQ-3000 Features & Options

### WORK HORSE OF THE AQUA PHASE LINE WITH

## DISINFECTANT APPLICATION SYSTEM

#### **Standard Features:**

- Double Wall Stainless Steel Chamber
- Stainless Steel Spray Arms
- Oscillating Spray Arms
- Wheelchair Wheel Rotation System
- Variable Wash Time Cycles (3, 5, or 7 Minute)
- Disinfectant Application System
- Wash and Rinse Cycles
- Automatic Chemical Injection
- LCD Status Display
- Front Door Window
- Heavy-Duty Caster Wheels
- Fill & Drain Hose Set
- 1.5 HP Motor/ Stainless Steel Pump
- 115V, 50/60Hz, 20 Amp AC Operation

#### **Optional Features:**

- Utility Equipment Rack
- Utility Wash Trays
- Air Dry Blower
- Lighted Wash Chamber
- 230V, 50/60Hz, Single Phase AC Operation
- Electronic Temperature Monitor



#### Total Outer Dimensions

• Width 33-1/2"

• Length 64"

• Height 62"

Effective Chamber Dimensions (with wheel rotation system)

• Width 26-1/2"

• Length 51-1/4"

• Height 36-3/4"

Like all Aqua Phase models, in order to operate the machine once the unit is located, all you need to do is hook up two hoses — one for the hot water fill and one for the drain — connect the Aqua Phase to an appropriate power source and begin saving time, improving efficiency and enhancing your infection control protocol.