



AQUA PHASE™

Washers & Sanitizers to automatically wash and apply disinfectants to wheelchairs and other hard-to-clean medical equipment.

AQ-2700

Features & Options

Stretched Version of the AQ-2000

Taller & Longer - see below

Standard Features:

- Double Wall Stainless Steel Chamber
- Stainless Steel Spray Arms
- Oscillating Spray Arms
- Variable Wash Time Cycles (3, 5, or 7 Minute)
- Wash and Rinse Cycles
- Automatic Chemical Injection
- LCD Status Display
- Heavy-Duty Caster Wheels
- Fill & Drain Hose Set
- 1 HP Motor / Stainless Steel Pump
- 115 Volt, 20 Amp AC Operation

Optional Features:

- Wheelchair Wheel Rotation System
- Disinfectant Application System
- Utility Equipment Rack
- Utility Wash Trays
- Front Door Window
- 1.5 HP Motor Upgrade / Stainless Steel Pump
- 230V, 50/60Hz, Single Phase AC Operation



Total Outer Dimensions

	<u>AQ-2000</u>	<u>AQ-2700</u>
• Width	33 1/2"	33 1/2"
• Length	64"	72"
• Height	62"	70 3/4"

Effective Chamber Dimensions

	<u>AQ-2000</u>	<u>AQ-2700</u>
• Width	26 1/2"	26 1/2"
• Length	51 1/4"	58 1/2"
• Height	36 3/4"	52"

Like all Aqua Phase models, in order to operate the machine once the unit is located, all you need to do is hook up two hoses — one for the hot water fill and one for the drain. Plug the Aqua Phase to an appropriate power source and begin saving time, improving efficiency and enhancing your infection control protocol.