



AQUA PHASE™

Washers & Sanitizers to automatically wash and apply disinfectants to wheelchairs and other hard-to-clean medical equipment.

AQ-3700

Features & Options

Stretched Version of the AQ-3000

Taller & Longer - see below

Standard Features:

- Double Wall Stainless Steel Chamber
- Stainless Steel Spray Arms
- Oscillating Spray Arms
- Variable Wash Time Cycles (3, 5, or 7 Minute)
- Disinfectant Application System
- Wash and Rinse Cycles
- Automatic Chemical Injection
- LCD Status Display
- Front Door Window
- Heavy-Duty Caster Wheels
- Fill & Drain Hose Set
- 1.5 HP Motor/ Stainless Steel Pump
- 115 Volt, 50/60Hz, 20 Amp AC Operation

Optional Features:

- Wheelchair Wheel Rotation System
- Utility Equipment Rack
- Utility Wash Trays
- Air Dry Blower
- Electronic Temperature Monitor
- 230V, 50/60Hz, Single Phase AC Operation



Total Outer Dimensions

| | <u>AQ-3000</u> | <u>AQ-3700</u> |
|----------|----------------|----------------|
| • Width | 33 1/2" | 33 1/2" |
| • Length | 64" | 72" |
| • Height | 62" | 70 3/4" |

Effective Chamber Dimensions

| | <u>AQ-3000</u> | <u>AQ-3700</u> |
|----------|----------------|----------------|
| • Width | 26 1/2" | 26 1/2" |
| • Length | 51 1/4" | 58 1/2" |
| • Height | 36 3/4" | 52" |

Like all Aqua Phase models, in order to operate the machine once the unit is located, all you need to do is hook up two hoses — one for the hot water fill and one for the drain. Plug the Aqua Phase to an appropriate power source and begin saving time, improving efficiency and enhancing your infection control